

Supplies for Candles & The Soap Kitchen Ltd Unit E Swinton Bridge Industrial Estate, White Lea Road,Swinton South Yorkshire S64 8BH, UK

Office Hours: Mon to Fri 08.00 to 16.00 customerservice@suppliesforcandles.co.uk

TEL: 01709 257151



Product Specification

Primesurf GMS SE40 Granules RSPO MB

INCI: Glyceryl Stearate

Description: Octadecanoic Acid, Monoester with 1,2,3-propanetriol

Stearic acid, Monoester with glycerol

CAS number: 31566-31-1 or 11099-07-3 or 85251-77-0 **EC number:** 250-705-4 or 234-325-6 or 286-490-9

Test parameter	Unit	Specification
Appearance		Small Beads/Granules
Acid Value	mg KOH/g	Max. 3
Iodine Value	g I2/100g	Max. 2
Free Glycerol	%	Max. 7
Total Monoglycerides	%	Min. 40
Melting point	°C	~58

Additional information

For additional information, please refer to the Material Safety Data Sheet.

Storage, packaging and handling

- Store in manufacturer's original, sealed container
- Store in a sheltered and well-ventilated warehouse away from moisture and sunlight
- Shelf life is 24 months from the date of manufacture

Disclaimer

Suggestion for the use of the application of the products and guide formulations are given just as information and without commitment. Such suggestions do not release customers from testing the products for their intended processes and purposes. The supplier does not assume any liability or risk involved in the use of its products as the conditions of use are beyond its control. The use of the products is solely responsible for compliance with all laws and regulations applying to the use of the products, including patents of third parties.